




Cannulated Screws, Staples





Cannulated Screws and Instrument Set 7,0 mm



T 500 Cannulated Screws and Instrument Set
T 111 Aluminium Case
T 501 Upper Tray
T 502 Lower Tray

Stainless Steel

Cat. No.:	Description	Quant. in Set
5-002	Cannulated screw \varnothing 7,0 mm / L 30 mm, 16 mm thread, standard	2
5-003	Cannulated screw \varnothing 7,0 mm / L 35 mm, 16 mm thread, standard	2
5-004	Cannulated screw \varnothing 7,0 mm / L 40 mm, 16 mm thread, standard	2
5-005	Cannulated screw \varnothing 7,0 mm / L 45 mm, 16 mm thread, standard	2
5-006	Cannulated screw \varnothing 7,0 mm / L 50 mm, 16 mm thread, standard	2
5-007	Cannulated screw \varnothing 7,0 mm / L 55 mm, 16 mm thread, standard	2
5-008	Cannulated screw \varnothing 7,0 mm / L 60 mm, 16 mm thread, standard	2
5-009	Cannulated screw \varnothing 7,0 mm / L 65 mm, 16 mm thread, standard	2
5-010	Cannulated screw \varnothing 7,0 mm / L 70 mm, 16 mm thread, standard	2
5-011	Cannulated screw \varnothing 7,0 mm / L 75 mm, 16 mm thread, standard	2
5-012	Cannulated screw \varnothing 7,0 mm / L 80 mm, 16 mm thread, standard	4
5-013	Cannulated screw \varnothing 7,0 mm / L 85 mm, 16 mm thread, standard	4
5-014	Cannulated screw \varnothing 7,0 mm / L 90 mm, 16 mm thread, standard	4
5-015	Cannulated screw \varnothing 7,0 mm / L 95 mm, 16 mm thread, standard	2
5-016	Cannulated screw \varnothing 7,0 mm / L 100 mm, 16 mm thread, standard	2
5-017	Cannulated screw \varnothing 7,0 mm / L 105 mm, 16 mm thread, standard	2
5-018	Cannulated screw \varnothing 7,0 mm / L 110 mm, 16 mm thread, standard	2
5-019	Cannulated screw \varnothing 7,0 mm / L 115 mm, 16 mm thread, standard	1
5-020	Cannulated screw \varnothing 7,0 mm / L 120 mm, 16 mm thread, standard	1
5-021	Cannulated screw \varnothing 7,0 mm / L 125 mm, 16 mm thread, standard	1
5-022	Cannulated screw \varnothing 7,0 mm / L 130 mm, 16 mm thread, standard	1

Cat. No.:	Description	Quant. in Set
5-031	Cannulated screw \varnothing 7,0 mm / L 45 mm, 32 mm thread, standard	2
5-032	Cannulated screw \varnothing 7,0 mm / L 50 mm, 32 mm thread, standard	2
5-033	Cannulated screw \varnothing 7,0 mm / L 55 mm, 32 mm thread, standard	2
5-034	Cannulated screw \varnothing 7,0 mm / L 60 mm, 32 mm thread, standard	2
5-035	Cannulated screw \varnothing 7,0 mm / L 65 mm, 32 mm thread, standard	2
5-036	Cannulated screw \varnothing 7,0 mm / L 70 mm, 32 mm thread, standard	2
5-037	Cannulated screw \varnothing 7,0 mm / L 75 mm, 32 mm thread, standard	2
5-038	Cannulated screw \varnothing 7,0 mm / L 80 mm, 32 mm thread, standard	2
5-039	Cannulated screw \varnothing 7,0 mm / L 85 mm, 32 mm thread, standard	2
5-040	Cannulated screw \varnothing 7,0 mm / L 90 mm, 32 mm thread, standard	2
5-041	Cannulated screw \varnothing 7,0 mm / L 95 mm, 32 mm thread, standard	1
5-042	Cannulated screw \varnothing 7,0 mm / L 100 mm, 32 mm thread, standard	1
5-043	Cannulated screw \varnothing 7,0 mm / L 105 mm, 32 mm thread, standard	1
5-044	Cannulated screw \varnothing 7,0 mm / L 110 mm, 32 mm thread, standard	1
5-045	Cannulated screw \varnothing 7,0 mm / L 115 mm, 32 mm thread, standard	1
5-046	Cannulated screw \varnothing 7,0 mm / L 120 mm, 32 mm thread, standard	1
5-047	Cannulated screw \varnothing 7,0 mm / L 125 mm, 32 mm thread, standard	1
5-048	Cannulated screw \varnothing 7,0 mm / L 130 mm, 32 mm thread, standard	1
5-522	Washer for screws \varnothing 6,5 / 7,0 mm	6

	1 Tray	
5-100	Drill Bit 5.0 mm \varnothing	1
5-114	Tap 7.0 mm \varnothing	1
5-118	Screw Driver	1
5-130	Tissue Protector Sleeve	1
5-140	Guide Wire 2.0 mm \varnothing 230 mm	10
4-184	Screw Holding Forceps	1
5-160	Direct Screws Length Gauge for \varnothing 7.0 mm	1
5-180	Cannulated Counter Sink	1
5-164	Canal Cleaner \varnothing 2 mm	1



Cannulated Screws



7 mm \varnothing with over flats SW 5.0 mm
16 mm length of thread, drill bit 5.0 mm
stainless steel

Cat. No.:	Screw Length
5-002	30 mm
5-003	35 mm
5-004	40 mm
5-005	45 mm
5-006	50 mm
5-007	55 mm
5-008	60 mm
5-009	65 mm
5-010	70 mm
5-011	75 mm
5-012	80 mm
5-013	85 mm
5-014	90 mm
5-015	95 mm
5-016	100 mm
5-017	105 mm
5-018	110 mm
5-019	115 mm
5-020	120 mm
5-021	125 mm
5-022	130 mm



7 mm \varnothing with over flats SW 5.0 mm
32 mm length of thread, drill bit 5.0 mm
stainless steel

Cat. No.:	Screw Length
5-031	45 mm
5-032	50 mm
5-033	55 mm
5-034	60 mm
5-035	65 mm
5-036	70 mm
5-037	75 mm
5-038	80 mm
5-039	85 mm
5-040	90 mm
5-041	95 mm
5-042	100 mm
5-043	105 mm
5-044	110 mm
5-045	115 mm
5-046	120 mm
5-047	125 mm
5-048	130 mm



7 mm \varnothing with over flats SW 5.0 mm
full threaded, drill bit 5.0 mm
stainless steel

Cat. No.:	Screw Length
5-270	45 mm
5-271	50 mm
5-272	55 mm
5-273	60 mm
5-274	65 mm
5-275	70 mm
5-276	75 mm
5-277	80 mm
5-278	85 mm
5-279	90 mm
5-280	95 mm
5-281	100 mm
5-282	105 mm
5-283	110 mm
5-284	115 mm
5-285	120 mm
5-286	125 mm
5-287	130 mm



Guide Wire with 10 mm thread

5-140 Guide Wire 2.0 mm \varnothing for Hollow Screws 7 mm \varnothing 230 mm


Instruments for 7 mm Cannulated Screws



5-100 Drill Bit \varnothing 5.0 mm 180 mm




5-180 Cannulated Counter Sink, Length 230 mm



5-114 Tap 7.0 mm \varnothing



5-118 5 mm Screw Driver for 7 mm screws



5-130 Tissue Protector Sleeve



5-160 Direct Screws Length Gauge for 7.0 mm \varnothing



5-162 Canal Brush 2.1 mm \varnothing



5-164 Canal Cleaner 2.0 mm \varnothing



5-178 Holding Clamp



Washers for Hollow Screws

5-132	19 mm \varnothing
5-134	16 mm \varnothing
5-522	13 mm \varnothing



5-131 Parallel Target Device for Guide Wires \varnothing 2.0 mm



T 510 Cannulated Screws and Instrument Set for 3.5 mm Screws



15-2162

L (mm)	B (mm)	H (mm)
310	190	65



Bottom perforated
Lid perforated

Set contents of implants

Cat. No.:	Description	Quant. in Set
5-300	Cannulated screw \varnothing 3,5 mm / L 10 mm, short thread, self-tapping	3
5-302	Cannulated screw \varnothing 3,5 mm / L 12 mm, short thread, self-tapping	3
5-304	Cannulated screw \varnothing 3,5 mm / L 14 mm, short thread, self-tapping	3
5-306	Cannulated screw \varnothing 3,5 mm / L 16 mm, short thread, self-tapping	3
5-308	Cannulated screw \varnothing 3,5 mm / L 18 mm, short thread, self-tapping	3
5-310	Cannulated screw \varnothing 3,5 mm / L 20 mm, short thread, self-tapping	3
5-312	Cannulated screw \varnothing 3,5 mm / L 22 mm, short thread, self-tapping	3
5-314	Cannulated screw \varnothing 3,5 mm / L 24 mm, short thread, self-tapping	3
5-316	Cannulated screw \varnothing 3,5 mm / L 26 mm, short thread, self-tapping	3
5-318	Cannulated screw \varnothing 3,5 mm / L 28 mm, short thread, self-tapping	3
5-320	Cannulated screw \varnothing 3,5 mm / L 30 mm, short thread, self-tapping	3
5-322	Cannulated screw \varnothing 3,5 mm / L 32 mm, short thread, self-tapping	3
5-324	Cannulated screw \varnothing 3,5 mm / L 34 mm, short thread, self-tapping	3
5-326	Cannulated screw \varnothing 3,5 mm / L 36 mm, short thread, self-tapping	3
5-328	Cannulated screw \varnothing 3,5 mm / L 38 mm, short thread, self-tapping	3
5-330	Cannulated screw \varnothing 3,5 mm / L 40 mm, short thread, self-tapping	3
5-332	Cannulated screw \varnothing 3,5 mm / L 42 mm, short thread, self-tapping	3
5-334	Cannulated screw \varnothing 3,5 mm / L 44 mm, short thread, self-tapping	3
5-336	Cannulated screw \varnothing 3,5 mm / L 46 mm, short thread, self-tapping	3
5-338	Cannulated screw \varnothing 3,5 mm / L 48 mm, short thread, self-tapping	3
5-340	Cannulated screw \varnothing 3,5 mm / L 50 mm, short thread, self-tapping	3

Cat. No.:	Description	Quant. in Set
5-342	Cannulated screw \varnothing 3,5 mm / L 10 mm, fully threaded, self-tapping	3
5-344	Cannulated screw \varnothing 3,5 mm / L 12 mm, fully threaded, self-tapping	3
5-346	Cannulated screw \varnothing 3,5 mm / L 14 mm, fully threaded, self-tapping	3
5-348	Cannulated screw \varnothing 3,5 mm / L 16 mm, fully threaded, self-tapping	3
5-350	Cannulated screw \varnothing 3,5 mm / L 18 mm, fully threaded, self-tapping	3
5-352	Cannulated screw \varnothing 3,5 mm / L 20 mm, fully threaded, self-tapping	3
5-354	Cannulated screw \varnothing 3,5 mm / L 22 mm, fully threaded, self-tapping	3
5-356	Cannulated screw \varnothing 3,5 mm / L 24 mm, fully threaded, self-tapping	3
5-358	Cannulated screw \varnothing 3,5 mm / L 26 mm, fully threaded, self-tapping	3
5-360	Cannulated screw \varnothing 3,5 mm / L 28 mm, fully threaded, self-tapping	3
5-362	Cannulated screw \varnothing 3,5 mm / L 30 mm, fully threaded, self-tapping	3
5-364	Cannulated screw \varnothing 3,5 mm / L 32 mm, fully threaded, self-tapping	3
5-366	Cannulated screw \varnothing 3,5 mm / L 34 mm, fully threaded, self-tapping	3
5-368	Cannulated screw \varnothing 3,5 mm / L 36 mm, fully threaded, self-tapping	3
5-370	Cannulated screw \varnothing 3,5 mm / L 38 mm, fully threaded, self-tapping	3
5-372	Cannulated screw \varnothing 3,5 mm / L 40 mm, fully threaded, self-tapping	3
5-374	Cannulated screw \varnothing 3,5 mm / L 42 mm, fully threaded, self-tapping	3
5-376	Cannulated screw \varnothing 3,5 mm / L 44 mm, fully threaded, self-tapping	3
5-378	Cannulated screw \varnothing 3,5 mm / L 46 mm, fully threaded, self-tapping	3
5-380	Cannulated screw \varnothing 3,5 mm / L 48 mm, fully threaded, self-tapping	3
5-382	Cannulated screw \varnothing 3,5 mm / L 50 mm, fully threaded, self-tapping	3
4-176	Washer 7,0 x 3,8 mm	5

9-334	Tap Sleeve 3,5 mm	
	Drill Sleeve 2,5 mm	1
9-167	Screw Driver Sleeve for 5-122	1
5-194	Tap for 3,5 mm cannulated Screws	1
5-188	Drill bit 3,5 mm / 1,4 mm cannulated	1
5-190	Drill bit 2,7 mm / 1,4 mm cannulated	1
9-093	Cannulated Tap Handle	1
5-126	Tap Sleeve 2,7 mm / Drill Sleeve 1,3 mm	1
5-184	Cannulated Counter sink	
	for 3,5 mm and 4,0 mm screws	1
5-122	Screw Driver f cannulated Screws	1
5-135	Guide Wire 1,2 x 150 mm	10

5-172	Direkt Screws length Gauge for 3,5 mm	1
5-174	Canal Brush \varnothing 1,35 mm	1
5-176	Canal cleaner \varnothing 1,2 mm	1
4-184	Screws forceps	1
9-116	Sharp Hook	1
5-178	Holding Clamps	1



Cannulated Screws \varnothing 3.5 mm



4-176
Washers
7,0 x 3,6 mm


Dia of thread	3.5 mm	3.5 mm
Dia of core	2,5 mm	2,5 mm
Dia of shaft	2,5 mm	
Dia of head	5,0 mm	5,0 mm
Pitch	1,25 mm	1,25 mm
Cannulation	1,3 mm	1,3 mm
Length		
10 mm	5-300	5-342
12 mm	5-302	5-344
14 mm	5-304	5-346
16 mm	5-306	5-348
18 mm	5-308	5-350
20 mm	5-310	5-352
22 mm	5-312	5-354
24 mm	5-314	5-356
26 mm	5-316	5-358
28 mm	5-318	5-360
30 mm	5-320	5-362
32 mm	5-322	5-364
34 mm	5-324	5-366
36 mm	5-326	5-368
38 mm	5-328	5-370
40 mm	5-330	5-372
42 mm	5-332	5-374
44 mm	5-334	5-376
46 mm	5-336	5-378
48 mm	5-338	5-380
50 mm	5-340	5-382



Guide Wire with 10 mm thread

5-135 Guide Wire 1.2 mm \varnothing for Hollow Screws 3.5 mm \varnothing 150 mm

Instruments for 3.5 mm Cannulated Screws



9-334 Tap Sleeve 3.5 mm
Drill Sleeve 2.5 mm



9-168 Screwdriver Sleeve for 5-122



5-122 Screw Driver for 3.5 mm and 4.0 mm cannulated screws 195 mm



5-194 Tap for 3.5 mm cannulated screws 160 mm



9-093 Can Tap Handle
Cannulated



5-188 Cannulated Drill Bit 160 mm length
Drill Bit ø 3.5 mm / 1.35 mm
5-190 Drill Bit ø 2.7 mm / 1.35 mm



5-126 Tap Sleeve 2.7 mm
Drill Sleeve 1.0 mm

Instruments for 3.5 mm Cannulated Screws



5-184 Cannulated Counter Sink for 3,5 mm and 4,0 mm screws



4-184 Screw Forceps



9-116 Sharp Hook
for checking of fracture ends and
for removing in grown tissue in screw heads




5-172 Direct Screws Length Gauge for 3.5 mm ø



5-174 Canal Brush 1.35 mm ø



5-176 Canal Cleaner 1.2 mm ø



5-178 Holding Clamp



T 313 Rack for ø 3.5 mm 1.25 pitch screws



T 520 Cannulated Screws and Instrument Set for 4.5 mm Screws



15-2164

L (mm)	B (mm)	H (mm)
310	190	130



Bottom perforated
Lid perforated

Set contents of implants

Cat. No.:	Description	Quant. in Set
5-200	Cannulated screw \varnothing 4,5 mm / L 20 mm, short thread, self-tapping	3
5-201	Cannulated screw \varnothing 4,5 mm / L 22 mm, short thread, self-tapping	3
5-202	Cannulated screw \varnothing 4,5 mm / L 24 mm, short thread, self-tapping	3
5-203	Cannulated screw \varnothing 4,5 mm / L 26 mm, short thread, self-tapping	3
5-204	Cannulated screw \varnothing 4,5 mm / L 28 mm, short thread, self-tapping	3
5-205	Cannulated screw \varnothing 4,5 mm / L 30 mm, short thread, self-tapping	3
5-206	Cannulated screw \varnothing 4,5 mm / L 32 mm, short thread, self-tapping	3
5-207	Cannulated screw \varnothing 4,5 mm / L 34 mm, short thread, self-tapping	3
5-208	Cannulated screw \varnothing 4,5 mm / L 36 mm, short thread, self-tapping	3
5-209	Cannulated screw \varnothing 4,5 mm / L 38 mm, short thread, self-tapping	3
5-210	Cannulated screw \varnothing 4,5 mm / L 40 mm, short thread, self-tapping	3
5-211	Cannulated screw \varnothing 4,5 mm / L 42 mm, short thread, self-tapping	3
5-212	Cannulated screw \varnothing 4,5 mm / L 44 mm, short thread, self-tapping	3
5-213	Cannulated screw \varnothing 4,5 mm / L 46 mm, short thread, self-tapping	3
5-214	Cannulated screw \varnothing 4,5 mm / L 48 mm, short thread, self-tapping	3
5-215	Cannulated screw \varnothing 4,5 mm / L 50 mm, short thread, self-tapping	3
5-216	Cannulated screw \varnothing 4,5 mm / L 52 mm, short thread, self-tapping	3
5-217	Cannulated screw \varnothing 4,5 mm / L 54 mm, short thread, self-tapping	3
5-218	Cannulated screw \varnothing 4,5 mm / L 56 mm, short thread, self-tapping	3
5-219	Cannulated screw \varnothing 4,5 mm / L 58 mm, short thread, self-tapping	3
5-220	Cannulated screw \varnothing 4,5 mm / L 60 mm, short thread, self-tapping	3


Cat. No.:	Description	Quant. in Set
5-230	Cannulated screw \varnothing 4,5 mm / L 20 mm, fully thread, self-tapping	3
5-231	Cannulated screw \varnothing 4,5 mm / L 22 mm, fully thread, self-tapping	3
5-232	Cannulated screw \varnothing 4,5 mm / L 24 mm, fully thread, self-tapping	3
5-233	Cannulated screw \varnothing 4,5 mm / L 26 mm, fully thread, self-tapping	3
5-234	Cannulated screw \varnothing 4,5 mm / L 28 mm, fully thread, self-tapping	3
5-235	Cannulated screw \varnothing 4,5 mm / L 30 mm, fully thread, self-tapping	3
5-236	Cannulated screw \varnothing 4,5 mm / L 32 mm, fully thread, self-tapping	3
5-237	Cannulated screw \varnothing 4,5 mm / L 34 mm, fully thread, self-tapping	3
5-238	Cannulated screw \varnothing 4,5 mm / L 36 mm, fully thread, self-tapping	3
5-239	Cannulated screw \varnothing 4,5 mm / L 38 mm, fully thread, self-tapping	3
5-240	Cannulated screw \varnothing 4,5 mm / L 40 mm, fully thread, self-tapping	3
5-241	Cannulated screw \varnothing 4,5 mm / L 42 mm, fully thread, self-tapping	3
5-242	Cannulated screw \varnothing 4,5 mm / L 44 mm, fully thread, self-tapping	3
5-243	Cannulated screw \varnothing 4,5 mm / L 46 mm, fully thread, self-tapping	3
5-244	Cannulated screw \varnothing 4,5 mm / L 48 mm, fully thread, self-tapping	3
5-245	Cannulated screw \varnothing 4,5 mm / L 50 mm, fully thread, self-tapping	3
5-246	Cannulated screw \varnothing 4,5 mm / L 52 mm, fully thread, self-tapping	3
5-247	Cannulated screw \varnothing 4,5 mm / L 54 mm, fully thread, self-tapping	3
5-248	Cannulated screw \varnothing 4,5 mm / L 56 mm, fully thread, self-tapping	3
5-249	Cannulated screw \varnothing 4,5 mm / L 58 mm, fully thread, self-tapping	3
5-250	Cannulated screw \varnothing 4,5 mm / L 60 mm, fully thread, self-tapping	3
3-520	Washer 11 x 4,8 mm	5

5-260	Guide Wire \varnothing 1,6 mm for hollow Screws \varnothing 4,5 mm, 230 mm length	10
9-330	Tap and Drill Sleeve 4,5 / 32 mm	1
5-120	Screw Driver cannulated for 4,5 mm Screws	1
5-192	Tap for 4,5 mm cannulated screws AO	1
5-187	Drill bit \varnothing 4,5 mm cannulated 1,8 mm	1
5-189	Drill bit \varnothing 3,2 mm cannulated 1,8 mm	1
9-093	Cannulated T-shaped handle with quick coupling	1
9-166	Screw Driver Sleeve f. 5-120	1

5-182	cannulated counter sink for 4,5 mm screws with quick coupling	1
4-184	screw forceps self holding	1
5-166	direkt screws length gauge	1
5-168	Canal brush \varnothing 1,75 mm	1
5-170	Canal cleaner 1,6 mm	1
5-178	holding clamp	1
9-334	Tap and drill sleeve 3,5 / 2,5 mm	1



Cannulated Screws \varnothing 4.5 mm



3-520
Washers
11,0 x 4,8 mm


Dia of thread	4.5 mm	4.5 mm
Dia of core	3,0 mm	3.0 mm
Dia of shaft	3,0 mm	
Dia of head	6,5 mm	6,5 mm
Cannulation	1,75 mm	1,75 mm
Hexagon	3,5	3,5 mm
Length		
20 mm	5-200	5-230
22 mm	5-201	5-231
24 mm	5-202	5-232
26 mm	5-203	5-233
28 mm	5-204	5-234
30 mm	5-205	5-235
32 mm	5-206	5-236
34 mm	5-207	5-237
36 mm	5-208	5-238
38 mm	5-209	5-239
40 mm	5-210	5-240
42 mm	5-211	5-241
44 mm	5-212	5-242
46 mm	5-213	5-243
48 mm	5-214	5-244
50 mm	5-215	5-245
52 mm	5-216	5-246
54 mm	5-217	5-247
56 mm	5-218	5-248
58 mm	5-219	5-249
60 mm	5-220	5-250



Guide Wire with 10 mm thread

5-260 Guide Wire 1.6 mm \varnothing for Hollow Screws \varnothing 4,5 mm, 230 mm

Instruments for 4.5 mm Cannulated Screws



9-334 Tap Sleeve 3.5 mm
Drill Sleeve 2.5 mm



9-166 Sleeve for 5-120 Driver



5-120 Screw Driver for 4.5 mm cannulated screws 225 mm




5-192 Tap for 4.5 mm cannulated screws 160 mm



9-093 Can Tap Handle
cannulated



5-187 Cannulated Drill Bit 160 mm length
Drill Bit \varnothing 4.5 mm / 1.8 mm
5-189 Drill Bit \varnothing 3.2 mm / 1.8 mm



9-330 Tap Sleeve 4.5 mm
Drill Sleeve 3.2 mm



Instruments for cannulated Screws



5-182 Cannulated Counter Sink for 4,5 mm screws



4-184 Screw Forceps



9-116 Sharp Hook
for checking of fracture ends and
for removing in grown tissue in screw heads




5-166 Direct Screws Length Gauge for 4.5 mm \varnothing



5-168 Canal Brush 1.75 mm \varnothing



5-170 Canal Cleaner 1.6 mm \varnothing



5-178 Holding Clamp



T 530 Cannulated Screws and Instrument Set for 4.0 mm Screws



15-2162

L (mm)	B (mm)	H (mm)
310	190	65



Bottom perforated
Lid perforated

Set contents of implants

Cat. No.:	Description	Quant. in Set
5-400	Cannulated screw \varnothing 4,0 mm / L 10 mm, short thread, self-tapping	3
5-401	Cannulated screw \varnothing 4,0 mm / L 12 mm, short thread, self-tapping	3
5-402	Cannulated screw \varnothing 4,0 mm / L 14 mm, short thread, self-tapping	3
5-403	Cannulated screw \varnothing 4,0 mm / L 16 mm, short thread, self-tapping	3
5-404	Cannulated screw \varnothing 4,0 mm / L 18 mm, short thread, self-tapping	3
5-405	Cannulated screw \varnothing 4,0 mm / L 20 mm, short thread, self-tapping	3
5-406	Cannulated screw \varnothing 4,0 mm / L 22 mm, short thread, self-tapping	3
5-407	Cannulated screw \varnothing 4,0 mm / L 24 mm, short thread, self-tapping	3
5-408	Cannulated screw \varnothing 4,0 mm / L 26 mm, short thread, self-tapping	3
5-409	Cannulated screw \varnothing 4,0 mm / L 28 mm, short thread, self-tapping	3
5-410	Cannulated screw \varnothing 4,0 mm / L 30 mm, short thread, self-tapping	3
5-411	Cannulated screw \varnothing 4,0 mm / L 32 mm, short thread, self-tapping	3
5-412	Cannulated screw \varnothing 4,0 mm / L 34 mm, short thread, self-tapping	3

Cat. No.:	Description	Quant. in Set
5-413	Cannulated screw \varnothing 4,0 mm / L 36 mm, short thread, self-tapping	3
5-414	Cannulated screw \varnothing 4,0 mm / L 38 mm, short thread, self-tapping	3
5-415	Cannulated screw \varnothing 4,0 mm / L 40 mm, short thread, self-tapping	3
5-416	Cannulated screw \varnothing 4,0 mm / L 42 mm, short thread, self-tapping	3
5-417	Cannulated screw \varnothing 4,0 mm / L 44 mm, short thread, self-tapping	3
5-418	Cannulated screw \varnothing 4,0 mm / L 46 mm, short thread, self-tapping	3
5-419	Cannulated screw \varnothing 4,0 mm / L 48 mm, short thread, self-tapping	3
5-420	Cannulated screw \varnothing 4,0 mm / L 50 mm, short thread, self-tapping	3
3-520	Washer 11 x 4,8 mm	5

9-334	Tap Sleeve 3,5 mm	
	Drill Sleeve 2,5 mm	1
9-168	Screw Driver Sleeve for 5-122	1
5-193	Tap for 4,0 mm cannulated screws quick coupling	1
5-188	Drill Bit 3,5 / 1,4 mm cannulated	1
5-190	Drill Bit 2,7 mm / 1,4 mm cannulated	1
9-093	Cannulated Tap handle	1
5-126	Tap Sleeve 2,7 mm Drill sleeve 1,3 mm	1
5-184	Cannulated Counter sink for 3,5 and 4,0 mm screws	1
4-184	screws forceps	1
9-116	sharp hook	1
5-172	Direkt screws length guage for 4,0 mm	1
5-174	Canal brush 1,35 mm	1
5-176	Canal cleaner 1,2 mm	1
5-178	Holding clamps	1
5-135	Guide Wire \varnothing 1,2 mm x 150 mm	10
5-122	Screw Driver f. cannulated screws	1



Cannulated Screws \varnothing 4.0 mm



3-520
Washers
11,0 x 4,8 mm


Dia of thread	4.0 mm
Pitch	1,75 mm
Dia of core	2,5 mm
Dia of shaft	2,7 mm
Dia of head	5,0 mm
Cannulation	1,35 mm
Hexagon	2,5
Length	
10 mm	5-400
12 mm	5-401
14 mm	5-402
16 mm	5-403
18 mm	5-404
20 mm	5-405
22 mm	5-406
24 mm	5-407
26 mm	5-408
28 mm	5-409
30 mm	5-410
32 mm	5-411
34 mm	5-412
36 mm	5-413
38 mm	5-414
40 mm	5-415
42 mm	5-416
44 mm	5-417
46 mm	5-418
48 mm	5-419
50 mm	5-420



Guide Wire with 10 mm thread

5-135 Guide Wire 1.2 mm \varnothing for Hollow Screws \varnothing 4,0 mm, 150 mm

Instruments for 4.0 mm Cannulated Screws



9-334 Tap Sleeve 3.5 mm
Drill Sleeve 2.5 mm



9-168 Screwdriver Sleeve for 5-122



5-122 Screw Driver for 3.5 / 4.0 mm cannulated screws 195 mm



5-193 Tap for 4.0 mm cannulated screws AO 160 mm



9-093 Can Tap Handle
cannulated



5-188 Cannulated Drill Bit 160 mm length
Drill Bit ø 3.5 mm / 1.35 mm
5-190 Drill Bit ø 2.7 mm / 1.35 mm



5-126 Tap Sleeve 2.7 mm
Drill Sleeve 1.0 mm

Instruments for 3.5 mm Cannulated Screws



5-184 Cannulated Counter Sink for 3,5 mm and 4,0 mm screws



4-184 Screw Forceps



9-116 Sharp Hook
for checking of fracture ends and
for removing in grown tissue in screw heads




5-172 Direct Screws Length Gauge for 3.5 mm ø



5-174 Canal Brush 1.35 mm ø



5-176 Canal Cleaner 1.2 mm ø



5-178 Holding Clamp



T 315 Rack for ø 4.0 mm screws (half threaded)



T 5400 Cannulated Screws and Instrument Set for 3.0 mm Screws TITAN



15-2162

L (mm)	B (mm)	H (mm)
310	190	65



Bottom perforated
Lid perforated

Set contents of implants TITAN

Cat. No.:	Description	Quant. in Set
5-50000	Cannulated screw \varnothing 3,0 mm / L 10 mm, short thread, self-tapping	3
5-50100	Cannulated screw \varnothing 3,0 mm / L 12 mm, short thread, self-tapping	3
5-50200	Cannulated screw \varnothing 3,0 mm / L 14 mm, short thread, self-tapping	3
5-50300	Cannulated screw \varnothing 3,0 mm / L 16 mm, short thread, self-tapping	3
5-50400	Cannulated screw \varnothing 3,0 mm / L 18 mm, short thread, self-tapping	3
5-50500	Cannulated screw \varnothing 3,0 mm / L 20 mm, short thread, self-tapping	3
5-50600	Cannulated screw \varnothing 3,0 mm / L 22 mm, short thread, self-tapping	3
5-50700	Cannulated screw \varnothing 3,0 mm / L 24 mm, short thread, self-tapping	3
5-50800	Cannulated screw \varnothing 3,0 mm / L 26 mm, short thread, self-tapping	3
5-50900	Cannulated screw \varnothing 3,0 mm / L 28 mm, short thread, self-tapping	3
5-51000	Cannulated screw \varnothing 3,0 mm / L 30 mm, short thread, self-tapping	3
5-51100	Cannulated screw \varnothing 3,0 mm / L 32 mm, short thread, self-tapping	3
5-51200	Cannulated screw \varnothing 3,0 mm / L 34 mm, short thread, self-tapping	3
5-51300	Cannulated screw \varnothing 3,0 mm / L 36 mm, short thread, self-tapping	3
5-51400	Cannulated screw \varnothing 3,0 mm / L 38 mm, short thread, self-tapping	2
5-51500	Cannulated screw \varnothing 3,0 mm / L 40 mm, short thread, self-tapping	2

5-186	Cannulated Countersink for 3,0 mm screws	1
9-093	Tap Handle	1
5-191	Cannulated Drill bit 2,0 / 1,2 mm	2
9-262	Drill Sleeve 2,0 mm	
	Guide Wire 1,0 mm	1
9-260	Drill Sleeve 5,5 / 4,3 mm	1
9-116	Sharp Hook	1
5-172	Direkt screws length gauge for 3,0 mm	1
5-174	Canal brush 1,25 mm	1
5-176	Canal cleaner 1,1 mm	1
5-178	Holding clamp	1
5-5400	Washer 6,5 / 3,3 mm	5
5-496	Kirschner Wire	10
4-184	Screw forceps	1
5-123	Screwdriver cannulated for 3.0 mm screws	1
5-125	Selfholding sleeve for Screwdriver 5-123	1



Cannulated Screws \varnothing 3.0 mm TITAN



5-5400

Washers
6,5 x 3,3 mm

Dia of thread	3.0 mm
Dia of core	2,0 mm
Dia of head	4,0 mm
Cannulation	1,2 mm
Length	
10 mm	5-50000
12 mm	5-50100
14 mm	5-50200
16 mm	5-50300
18 mm	5-50400
20 mm	5-50500
22 mm	5-50600
24 mm	5-50700
26 mm	5-50800
28 mm	5-50900
30 mm	5-51000
32 mm	5-51100
34 mm	5-51200
36 mm	5-51300
38 mm	5-51400
40 mm	5-51500



Guide Wire with 5 mm thread

5-496 Guide Wire 1.0 mm \varnothing for Hollow Screws \varnothing 3,0 mm, 150 mm



Instruments for 3.0 mm Cannulated Screws



5-186 Cannulated Counter Sink for 3,0 mm screws



4-184 Screw Forceps



9-116 Sharp Hook
for checking of fracture ends and
for removing in grown tissue in screw heads



5-191 Cannulated Drill Bit
Drill Bit \varnothing 2.0 mm / 1.2 mm



9-262 Drill Sleeve Drill Bit 2.0 mm
and Guide wire 1.0 mm




9-260 Drill Sleeve 5.5 mm / 4,3 mm

Instruments for 3.0 mm Cannulated Screws



5-172 Direct Screws Length Gauge for 3.5 mm \varnothing



5-174 Canal Brush 1.35 mm \varnothing



5-176 Canal Cleaner 1.2 mm \varnothing



5-178 Holding Clamp



T 313 Rack for \varnothing 3.0 mm screws (half threaded)



9-093 Tap Handle for AO shaft
cannulated



5-123 Screwdriver for 3.0 mm screws
cannulated



5-125 Self holding Sleeve for Screwdriver 5-123

Epiphyseal Staples Stainless Steel



5-631
For conveytry



- 5-600** Barbed Staples 8.7 x 16 x 2 mm
- 5-602** Barbed Staples 8.7 x 19 x 2 mm
- 5-604** Barbed Staples 8.7 x 22 x 2 mm
- 5-606** Barbed Staples 8.7 x 25 x 2 mm



5-610
24 x 33 x 4.8 mm



5-612
24 x 33 x 9.5 mm



5-614
24 x 33 x 14 mm



	Diameter	Width	Length
5-620	2.5 mm	22 mm	19 mm
5-622	2.5 mm	16 mm	19 mm
5-624	1.5 mm	22 mm	19 mm
5-626	1.5 mm	16 mm	19 mm



5-630 Staples Extractor



5-632 Staples Driver



5-634 Staples Inserter for U-Staples



Staples



Chron Cobalt Molybden
5-650 6 mm



Chron Cobalt Molybden
5-652 8 mm



Chron Cobalt Molybden
5-654 11 mm



Chron Cobalt Molybden
5-656 14 mm

Titan
5-65000 6 mm

Titan
5-65200 8 mm

Titan
5-65400 11 mm


Titan
5-65600 14 mm



1-1920 Schlitzhammer mit Schlitzweite 14 mm
Slotted Hammer - width of the slot 14 mm



Chrom Cobalt Molybden			
Cat. No.:	Width	Length	Height
5-730	24 mm	40 mm	5 mm
5-732	24 mm	40 mm	10 mm
5-734	24 mm	40 mm	15 mm




Cat. No.: 5-670 Impactor-Extractor for Staples 5-650 - 5-656 and 5-730 - 5-734


5-670 Impactor




Chrome Cobalt Molybden Osteotomy Staples



Cat. No.:	Width A	Length B	Height C
5-710	24 mm	30 mm	5 mm
5-712	24 mm	30 mm	10 mm
5-714	24 mm	30 mm	15 mm



Cat. No.:	Width	Length
5-705	16 mm	22 mm
5-706	22 mm	22 mm
5-707	28 mm	22 mm





Interferenzschrauben kannuliert
Interference screws cannulated



Material: Titan ISO 5832 / 03

Total Length / mm	20	25	30
Flat Head			
ø 7 mm	5-900	5-902	5-904
ø 8 mm	5-910	5-912	5-914
ø 9 mm	5-920	5-922	5-924



Total Length / mm	20	25	30
Round Head			
ø 7 mm	5-927	5-928	5-929
ø 8 mm	5-930	5-931	5-932
ø 9 mm	5-933	5-934	5-935

Cannulated



Cannulated Drill

5-950	ø 6.0 mm
5-952	ø 7.0 mm
5-954	ø 8.0 mm
5-956	ø 9.0 mm
5-958	ø 10.0 mm
5-960	ø 11.0 mm
5-962	ø 12.0 mm



5-120 Screw driver cannulated for interference screws

Interferenzschrauben, Schutzhülse
Interference screws, Protective Sleeves



7.0 mm	8.0 mm	9.0 mm
5-980	5-982	5-984

Tendon Button

4 Loch / Holes
4 x 12 mm



18-001



Suture Disk, \varnothing 14 mm
Reintitan
Pure titanium
18-030



18-040

Passierdraht \varnothing 2,4 mm mit Öse und Spiralspitze Länge: 28 cm
Passing Pin \varnothing 2.4 mm with Eye and Drill Tip Length: 28 cm



Suture Disk, \varnothing 14 mm
mit Rotationssperre
Reintitan
with rotation fixation
Pure titanium
18-032

